


A Cross-Modal ECG–SQUID MCG Representation Learning Framework

Rafsan Jany
Korean Convergence Medical Science

What is SQUID Magnetocardiography (MCG)

MCG = Magnetic ECG measured using SQUID devices

- SQUID = Superconducting Quantum Interference Device
- Most sensitive magnetic sensor in the world
- Detects magnetic fields produced by the heart
- Unlike ECG, magnetic fields **aren't distorted** by body tissues
- Can reveal cardiac abnormalities that ECG may miss
- Used today mainly in research centers (rare & expensive)



- Fenici et al., *European Heart Journal*, 2020
- Park et al., *Sensors*, 2020
- Koch et al., *Biomedical Engineering Online*, 2011


Why MCG matters (Clinical motivation)

MCG can detect:

- Ischemia earlier than ECG
- Conduction abnormalities
- Fetal arrhythmias (when ECG is unreliable)
- Subtle repolarization changes
- Spatial current distribution of the heart

Problem:

- Almost no hospitals use MCG – devices are rare, expensive, and datasets are tiny.
- This makes AI or deep-learning on MCG nearly impossible.



Gap 1 – MCG Technology Is Underutilized in Clinical Practice

- SQUID-MCG can detect ischemia, fetal arrhythmia, and subtle conduction delays
- But MCG is rare, expensive(*cost \$1M+*), and only available in research centers

Gap 2 – ECG Alone Cannot Capture Everything

- ECG signals get distorted by tissue, muscle, fat
- MCG is cleaner, but too rare
- Clinicians lack a unified view combining both modalities

Gap 3 – There Are Almost No MCG Datasets for AI

- Only one open-access dataset contains real SQUID-MCG: Koch et al. (2011)
- It is extremely small (127 beats)
- Makes modern AI training impossible without new methods

Gap 4 – No way to help hospitals without SQUID machines

- If we could use ECG to approximate MCG
- Many hospitals could benefit from MCG insights
- But no such system exists today

Probable Solution

We need a way to learn MCG representations without needing thousands of MCG recordings.

ECG and MCG measure the same heart activity, but each sees a DIFFERENT part of the truth.

- ECG beat → vector (electrical potential ON the skin)
- MCG beat → vector (magnetic field generated by the SAME currents INSIDE the heart)
- Both vectors are close if they come from the same heartbeat.

Simple Word:

ECG and MCG are like two cameras taking a picture of the same object from different angles.

- Both capture the heartbeat.
- Both look at the same depolarization/repolarization.
- Both detect conduction abnormalities.

BUT they see different distortions and different hidden details.

ECG → foundational cardiac patterns

MCG → rare, high-value complementary info

Core Research Idea

This study proposed: A model connects electrical and magnetic views of the same heartbeat to provide a more complete cardiac understanding.

The input to the trained model can be either an ECG alone, an MCG alone, or both – because the model learns a shared representation for both modalities. A shared latent heartbeat representation for clinical use.

PTB Diagnostic ECG Dataset

Source: PhysioNet / PTB Diagnostic ECG Database

Signals: 12-lead ECG

Sampling Rate: 1000 Hz

Subjects: 290+

Conditions: Myocardial infarction, arrhythmia, block, hypertrophy, healthy controls

Size: Thousands of heartbeats (large dataset)

PTB gives us a strong, robust ECG understanding before we try to align ECG with MCG.

Koch ECG-MCG (SQUID) Paired Dataset

Source: Koch et al., Biomedical Engineering Online (2011)

Modalities:

BSPM (Body Surface Potential Mapping) → electrical activity

- 32 channels
- 1000 Hz

MCG (Magnetocardiography) using SQUID

- 100 channels
- 1000 Hz
- Femto-Tesla units (fT)

MRI torso scans

Methodology Overview


Stage 1 – Learn ECG features from a large dataset (PTB Diagnostic ECG)

- Self-supervised learning (no labels needed)
- Model learns patterns of electrical activity

Stage 2 – Align ECG with MCG using paired data (Koch dataset)

- For each heartbeat: ECG signal \leftrightarrow MCG signal
- Use contrastive learning to match them in a shared space

Stage 3 – Align ECG with MCG using paired data (Koch dataset)



In this project we propose to:

Self-supervised learning (no labels needed)

- Build the first cross-modal representation between ECG and MCG
- Overcome MCG data scarcity using contrastive learning
- Explore a new direction for cardiac AI research
- Enable practical future applications like virtual MCG, enhanced diagnosis, and multimodal fusion

Status:

- Dataset (ready)
- Preprocessing complete (initially done)
- Model design complete (On progress: experimented with 1d CNN)
- Training (Next Step)
- Evaluation (Next Step)

Thank You